

SCPline

Make single-cell analysis available to every investigator

SCPline: An interactive shiny framework for the study of single-cell proteomics data preprocessing

SCPline is an integrated one-stop single-cell proteomics data preprocessing analysis tool developed by Guohua Wang Laboratory, aiming to provide efficient data preprocessing analysis for laboratory researchers.

Moudle2:scMultimics workflow

Multi-omics

Overview

Upload your data

Normalization & Feature Selection

PCA

UMAP

Clustering

Download Results

Load example data

Load Example and run

Upload your input data

SingleCellExperiment Object (Accepted Format: .rds)

Browse... No file selected

Analysis method: Seurat

Min nFeature:

200

Min nCells:

3

Project Name:

Project1

Load Data and run

Reset



You can try loading sample files or you can load your own data

Min nFeature: Sets the minimum feature filtering standard

Min nCells: Set the minimum number of cells for filtering

Project Name: Enter a project name for management and identification analysis.

Click the "Load Data and run" button to start data loading and analysis, or use the "Reset" button to reset the input.



Multi-omics

Overview

Upload your data

Normalization & Feature Selection

PCA

UMAP

Clustering

Download Results

Data Preprocessing

Number of variable features:

2000


Variable feature selection method:

vst

Identify highly variable features

Number of variable features: Sets the number of highly variable features to be selected, for example 2000.

Variable feature selection method: Selects a method for feature selection, for example, vst.



Multi-omics

Overview

Upload your data

Normalization & Feature
Selection

PCA

UMAP


Clustering

Download Results



Principal Component Analysis (PCA)

Run PCA

Click on 'Run PCA' under the 'PCA' tab to visualize PCA plot.



Moudle2:scMultimics workflow



Next comes UMAP and clustering, where the cluster analysis is complete

Multi-omics

Overview

Upload your data

Normalization & Feature Selection

PCA


UMAP

Clustering

Download Results

Clustering

Resolution:



Run Clustering

Use the slider to adjust the clustering resolution to control the fineness of the clustering. Each time the parameter is selected, the corresponding clustering effect of the parameter is found in the drop-down box below

RNA Clustering

Group by for RNA:

orig.ident

The drop-down box contains the metadata attributes of the uploaded data, and the attribute content of each data can be intelligently identified

RNA Clustering



Multi-omics

Overview

Upload your data

Normalization & Feature Selection

PCA

UMAP

Clustering

Download Results

Visualization

Download RDS file (SingleCellExperiment Object)

Finally, click the download button to get the data after the pre-processing process

Metadata Information

orig.ident	nCount_ADT	nFeature_ADT	nCount_RNA	nFeature_RNA	sample	timepoint	UMAP_1	UMAP_2	percent.mito
Dx	4765.00	19	5002.00	1683	Pt-12-Dx	Dx	-2.64	-2.91	0.00
Dx	2031.00	19	7356.00	2353	Pt-12-Dx	Dx	3.13	-4.47	0.00
Dx	1799.00	19	5419.00	1754	Pt-12-Dx	Dx	2.97	-4.20	0.00
Dx	2715.00	19	6052.00	2068	Pt-12-Dx	Dx	4.47	-4.59	0.00
Dx	2734.00	19	8756.00	2799	Pt-12-Dx	Dx	1.50	-4.74	0.00
Dx	2278.00	19	6169.00	1909	Pt-12-Dx	Dx	-0.24	-4.55	0.00
Dx	4573.00	19	4068.00	1477	Pt-12-Dx	Dx	-2.38	-5.06	0.00
Dx	1585.00	19	9062.00	2507	Pt-12-Dx	Dx	3.23	-4.68	0.00
Dx	3683.00	19	10217.00	2322	Pt-12-Dx	Dx	-1.07	-3.54	0.00
Dx	900.00	19	8936.00	2275	Pt-12-Dx	Dx	-12.77	0.24	0.00